

TRIAD Treatment chair for Triage, Fast Track, and Vertical Mode

Reduces handling time while improving comfort and safety.

The chair features comfortable and durable cushions designed to reduce pressure points while maintaining lumbar support. The seat is also designed to help the patient maintain a good position during treatment. The cushions are fully moulded and seamless with injected viscose foam.

Cleaning: The vinyl upholstery is tested for common cleaning products used in hospitals and clinics. It is resistant to Oxivir and Virox; imperative when infection control is a priority.



Rechargeable battery pack.



Quickly reposition the patient in Triage or adjust the patient in Vertical Mode

Standard colors:

- Sky Blue
- Gunmetal

Special order:

- Cobalt
- Wedgewood
- Desert Rose
- Ash



Emergency foot controller to position the patient in Trendelenburg while keeping hands free.



Trendelenburg -10° is effortless and requires only 8.8 sec.



Retractable tray table positioned to each armrest.



Footrest cushion on the adjustable foot rest for smaller people.



The arms swivel up, away from the patient, or toward the patient



Mini exercise bike which can be mounted on the footrest to exercise during treatment.

TRIAD

Treatment chair for Triage, Fast Track, and Vertical Mode

TRIAD configuration

Included features:

- » Fixed transport handles
- » Electrically adjustable footrest
- » 5" Twin wheels with central lock
- » Hand controller
- » Pivoting armrests, synchronized with backrest
- » Electric Trendelenburg 8.8 sec.
- » Seamless moulded cushions
- » Moulded viscose foam pillow

TRIAD Comfort Options

Included features:

- » TRIAD base configuration
- » 17" x 13" Swivel table that can be installed on right or left of patient
- » Emergency foot controller for safe and hands-free Trendelenburg in 8.8 seconds
- » Rechargeable battery pack

Additional options available

- » SEP-01: Closure plates under armrests
- » SHELF-TAB: Retractable tray table 15.2 x 30.5 cm / 6" x 12" positioned on each armrest
- » FC-2: Rear central lock pedals only
- » FRT-2: Credit for removing electric footrest
- » CU6-FR-TRIAD: Footrest cushion for short people
- » MEBT: Mini exercise bike which can be mounted on the footrest to exercise during treatment

Specifications

Adjustments

- Height: 20" to 28"
- Backrest: 0° to 80°
- Footrest: 6"
- Legs: 0° to -50°
- Trendelenburg: -10°

Maximum weight capacity

- 205 kg / 450 lb using hospital-grade electrical components
- 120V - 3.9 AMP - 60 Hz

Dimensions

- Standard: 61 x 195.6 cm / 24" x 77"
- Armrests: 48.9 x 18.4 cm / 19.25" x 7.25"

Cushions

- Molded for comfort and durability
- Seamless to help prevent infection
- Vinyl upholstery resistant to Oxivir and Virox

Controls

- Hand controller

Wheels

- 12.7 cm / 5" Twin wheels with central lock



Comfortable and Flexible

TSK PRODUCTS

sales@tskproducts.com | www.tskproducts.com



Entirely designed and manufactured in Canada, locally manufactured offers you easy access to spare parts and customer service in English and French.

 **IBIOM**

BY  **umanoMedical**

L'Islet QC G0R 2C0 CANADA